

**QMC Instruments Ltd**

# **Model TK1810 Cryostat**



## **Operating Manual**

Date

## CONTENTS

**First – A Word of Warning**

**Introduction**

**Packing List**

**1. Unpacking and Preparing the System for Operation**

Initial inspection

Cryostat inspection

**2. Evacuating the Cryostat**

**3. Liquid Nitrogen Pre-Cool**

A word about your vacuum pump

The need to pre-cool with liquid nitrogen

Safety valves

The pre-cool period

Removing the liquid nitrogen from the reservoir

**4. Liquid Helium Transfer**

The liquid helium transfer tube

Helium gas recovery

Keeping the cryostat cold

**5 Cryostat Cryogenic Performance**

**Appendix A. Cryostat Test Log Sheet**

## First - A Word of Warning

### Using Cryogenics

Cryogenic liquids are potentially dangerous. If you are not already familiar with the standard procedures appropriate for the use of liquid nitrogen and liquid helium, please seek advice before proceeding.

Operating this equipment involves the use of vacuum and cryogenic liquids. Please read this manual carefully before you operate the cryostat. Although this is not a safety instruction manual, the text describes our own procedures and this may help to avoid accidents.

The photo below shows part of a damaged cryostat. We do not want this to happen to you. Please ensure that all personnel involved in the use of the cryostat are fully accustomed with the techniques involved.



## Introduction

This is a QMC Instruments Ltd. type TK1810 liquid helium cryostat built to our specification by our sister company Thomas Keating Ltd.

Item	Serial Number
TK1810 cryostat	XXXXX-X / XXXX

## Packing List

The following items are included in this shipment. Please check the contents against this list and contact QMC Instruments as soon as possible if you suspect that any items are damaged or missing.

### **Type TK1810 Thomas Keating Ltd cryostat**

- o Cryostat fitted with:
  - Over-pressure relief valve fitted to the cryostat top plate for safety
  - Non-return valve
- o Cryostat central neck safety baffle in place which includes:
  - Over-pressure relief valve
- o Outer vacuum case base plate with O-ring (fitted)
- o Base plate for the gas cooled radiation shield (fitted)
- o Spares kit which includes:
  - Set of screws
  - O-rings
  - M3 and M4 Allen keys
- o Cryostat operating manual

## 1. Unpacking and Preparing the System for Operation

Photographs included in this manual are general photos that may not be specific to your particular cryostat.

### Initial inspection


Please inspect the box in which the goods were shipped, and the contents, for any obvious sign that damage has occurred in transit. If you think that the package has been damaged in some way, please contact us before proceeding further. Your equipment is guaranteed for two years against failure resulting from effects beyond your control, and we will be happy to make any repairs at no cost to you during this time.

There is a spares kit accompanying the cryostat that contains O-rings, bolts, screws etc.

### Cryostat inspection


Refer to **Photos 1.1, 1.2 & 1.3**

To allow access to the bottom-plate invert the cryostat so that it rests on the stainless lifting ring. To avoid marking the ring, stand the cryostat on something to protect it such as soft tissue, cloth or bubble wrap.



**Photo 1.1.** View inside a TK1810 cryostat showing a detector mounted in position, and the three aluminium transit posts

The TK1810 cryostat has a radiation shield base-plate and a room temperature outer vacuum casing (OVC) base-plate. The radiation shield base plate is located using the set of M3 screws provided. The black OVC base-plate is located using the M4 socket-headed screws provided. It is important to check that the O-ring is in place, that it is clean, greased and that its seating is free of marks and scratches. The screws locating the OVC base-plate should not be over-tightened because this can distort the O-ring and may cause vacuum leaks. If the screws are equally tightened, it is normal for a small gap to show between the lip of the OVC base-plate and the bottom of the cryostat casing.



**Photo 1.2.** Cryostat base-plates



**Photo 1.3.** Radiation shield base-plate in position


## 2. Evacuating the Cryostat

Please refer to **photo 2.1**. Before cooling the cryostat, the vacuum chamber must be evacuated by connecting a suitable pump to the evacuation port located on the top-plate. The pump should be capable of reducing the pressure in the cryostat to below  $10^{-1}$ mbar. This can with time be achieved by using a rotary pump only, but for optimum cryogenic performance of the cryostat it is better to use a diffusion or turbo-molecular pump to reduce the pressure still further.

The pumping system should ideally have a pressure gauge measuring the pressure as close to the cryostat as possible. Always check the quality of the pump system and pumping line prior to opening the cryostat vacuum valve.

The vacuum valve should be opened very slowly when the pressure in the cryostat is at or close to atmospheric pressure. This prevents rapid pressure changes that risk damage to the delicate components inside the cryostat.

Typically, the cryostat could be ready for pre-cooling (refer to Section 3) after pumping for thirty minutes using a two stage pumping station.



**Photo 2.1.** Cryostat top-plate fittings



## 3. Liquid Nitrogen Pre-Cool

**IMPORTANT:** Refer to the warning at the front of the manual before proceeding with cryogenic cooling of this cryostat.

### A word about your vacuum pump

The pressure in the cryostat should drop rapidly when filling with liquid nitrogen because some of the gas, mainly oxygen, begins to cryopump (condense onto the cold surfaces). The cryostat can remain attached to the pump during the pre-cool period if the pump you are using is an oil diffusion or turbomolecular type pump with a base pressure lower than  $10^{-6}$  mbar. If you are only using a rotary pump, then the pressure in the cryostat will be lower during the pre-cool period than the pump is capable of generating, and the pump must therefore be detached immediately prior to cooling.

### The need to pre-cool with liquid nitrogen

When a satisfactory pressure has been reached in the cryostat vacuum chamber, it is necessary to pre-cool the cryostat with liquid nitrogen before cooling with liquid helium. This will reduce the amount of liquid helium used.


The cryostat neck baffle assembly, **photo 3.1**, should be unscrewed and removed from the central port to enable liquid nitrogen to be poured into the reservoir, **photos 2.1 and 3.2**.

For preference, transfer the liquid nitrogen directly from a pressurised liquid nitrogen storage Dewar which should take less than 10 minutes to complete. Alternatively, pour the liquid nitrogen using a bucket and a funnel, as depicted in **photo 3.2**. In this case, the funnel must be attached to a pipe which extends down into the neck and well into the reservoir itself. For a TK1810 cryostat a length of at least 180mm is needed. The pipe diameter should be about 6mm (1/4 inch) to allow both reasonable throughput and space outside of the pipe for boiling nitrogen gas to escape.

### Safety valves

The top-plate fittings are shown in **photo 2.1**. The reservoir access port should always be fitted with the non-return valve to stop the condensation of moisture within the neck. This moisture could freeze and block the neck of the cryostat which in turn could lead to failure and damage.

The cryostat neck baffle is shown in **photo 3.2**. The baffle incorporates an overpressure release valve. Should an ice blockage form in the central neck of the cryostat, gas will be unable to escape through the non-return valve. Such an event will cause the overpressure relief valve, located at the top of the baffle to open, thereby releasing pressure from the reservoir.



**Photo 3.1.** Cryostat neck baffle



**Photo 3.2.** Using a funnel to fill the cryostat with liquid nitrogen

## **The pre-cool period**

This cryostat is designed for use within a typical working day. It should be ready for use within 2 hours of connecting to a pumping station and should remain cold for the remainder of the working day. Note however that the length of pre-cool period will determine the efficiency of use of liquid helium. A minimum of 30 minutes pre-cooling with liquid nitrogen is recommended.

## **Removing the liquid nitrogen from the reservoir**

When the pre-cool period is complete the liquid nitrogen in the reservoir should be removed. This is done by first removing the neck baffle and then carefully inverting the cryostat and pouring the nitrogen into a suitable cryogen-safe container

It is important that all the liquid nitrogen is removed from the central reservoir before the liquid helium transfer is started. Any liquid nitrogen remaining in the central reservoir will be frozen by the liquid helium. Nitrogen ice forms an effective insulating layer which will prevent the cold plate reaching its intended operating temperature. A large amount of expensive liquid helium will also be wasted in creating a small amount of very cold nitrogen ice!

## 4. Liquid Helium Transfer

When you are certain that all the liquid nitrogen has been removed from the reservoir the cryostat can be filled with liquid helium. The cryostat should be arranged such that the transfer tube reaches the bottom of the cryostat and the storage Dewar simultaneously.

It is wasteful to transfer liquid helium too quickly. A rubber bladder can be used to control the pressure driving the transfer, and the rate of filling can be judged from the size of the plume of exhaust helium gas rising from the cryostat in the case of an open liquid helium transfer.

### The liquid helium transfer tube

It is important that the liquid helium transfer tube used is designed to suit both the cryostat and the liquid helium storage Dewar. The delivery end of the transfer tube should have a fully evacuated section with diameter approximately 6mm ( $\frac{1}{4}$  inch) and length at least 200mm. It should therefore permit liquid helium to be delivered efficiently into the central reservoir while at the same time leave space for spent helium gas to escape without a build-up of pressure within the cryostat.

QMC Instruments Ltd. can arrange to supply a suitable liquid helium transfer tube for your cryostat. We offer a rigid transfer tube, product code QTT/R, and a flexible transfer tube, product code QTT/F, with a reach in excess of 1000mm. Please contact us, or your supplier, if you have any questions regarding the suitability of your equipment.

**Photo 4.1** depicts a liquid helium transfer in progress. **Photo 4.2** shows a typical boil-off plume in the phase when the cryostat is cooling between 77K and 4.2K. **Photo 4.3** shows the larger, cloudier and more erratic plume, which results when the liquid helium reservoir is full. At this stage the transfer should be terminated. It should take about fifteen minutes for a TK1810 cryostat to cool down from 77K to 4.2K and to fill with liquid helium; and the whole process should consume about three litres of liquid helium.

### Helium gas recovery

Here in Cardiff we have no facilities for recovering spent helium gas, hence all the liquid helium transfers undertaken in our laboratories are “open” in the manner shown in the photos. However some installations offer recovery facilities whereby a helium return line is attached to the exhaust port of the cryostat. Use the black anodized aluminium tightening ring and O-ring from the central neck fitting to make a seal around the liquid helium transfer tube. Under such circumstances, a coarse flowmeter could be inserted in the return line to indicate flow rate from the transfer. Usually a steady flow-rate is indicated during the cool and fill phases of the transfer. When the reservoir is full however, the flow rate becomes erratic, and the transfer should be terminated. When the transfer is complete the transfer tube should be removed carefully but swiftly and the safety valves fitted without delay.



**Photo 4.1.** Liquid helium transfer



**Photo 4.2.** Helium gas exhaust during fill



**Photo 4.3.** Helium plume when complete

### **Keeping the cryostat cold**

It is important to keep all the neck fittings and safety valves in place whenever the cryostat is cold. If these are removed for liquid helium transfer, they should be removed only at the last moment when all other preparations have been made. They should be replaced as soon as the transfer tube is removed.

The first fill liquid helium hold-time may be shorter than that for subsequent helium fills. This is because the initial liquid helium boil off rate may be high if significant further cooling takes place when the transfer is complete.

The cryostat can be kept continuously cold by repeatedly replenishing the liquid helium. The subsequent fill hold time for the liquid helium is shown in **Table 5.1** in **Section 5**. While this TK1810 cryostat is principally designed to be used during the working day and to warm up overnight, it is possible to keep it cold overnight for continued use the following day. To do this, the cryostat should have been cooled during the day and another (subsequent) liquid helium transfer carried out late in the day to fill up the reservoir with liquid helium. It might be necessary to carry out this subsequent fill during the evening and to arrive early the next day to top up the cryostat again (refer to **Table 5.1, Section 5**), but it will allow for continued operation of the cryostat where time is at a premium.

When transferring liquid helium into a cryostat that already contains liquid helium, the transfer tube should be fully cooled before inserted into the cryostat neck. This prevents the warm transfer tube and warm helium gas from boiling away excessive amounts of the liquid helium already in the

cryostat. In this case the transfer tube is inserted into the storage Dewar and the pressure control bladder inflated slightly to pass gas through the tube to cool it. When the transfer tube has cooled, dense, milky helium gas emerges from the delivery end, **photo 4.4**, and the transfer tube can then be manoeuvred carefully to the cryostat and lowered into the central neck. The refill can then proceed in the way described above.

When the liquid helium transfer is complete we advise that you wait about fifteen minutes before using the cryostats to allow any components on the cold plate to reach thermal equilibrium.



**Photo 4.4.** Liquid helium emerging from a cold tube

## 5. Cryostat Cryogenic Performance

The liquid helium hold-time of the cryostat is measured in QMC Instruments Ltd. tests and tabulated in Table 5.1. The liquid helium boil-off is measured until the radiation shield within the cryostat reaches thermal equilibrium. When equilibrium is reached the base boil-off is measured and used to determine the liquid helium hold-time of the cryostat. Note that a first fill will not last as long as a subsequent fill hold-time due to the high initial boil-off when the cryostat is cooling from liquid nitrogen temperature.

In order to achieve these figures it is important that the operating instructions laid out in this manual are followed, and that care is taken to ensure that the cryostat is completely full before the liquid helium transfer is terminated.

Reservoir capacity / litres	1.0
Base helium boil-off / litres of gas per min at STP	1.1
Liquid helium first fill hold-time / hours	$11 \pm 2$
Liquid helium subsequent fill hold-time / hours	$17 \pm 2$

**Table 5.1.** System cryogenic performance of the TK1810 cryostat



## Contract details and guarantee

This equipment is guaranteed for a period of two years from the date of delivery against failure caused by defective materials or workmanship. Defective parts will be repaired or replaced on return to the final supplier at no cost, provided that failure is not due to misuse or mishandling after delivery. QMC Instruments Limited will assume no liability for loss of life or damage to property arising from the use or misuse of its products.

Purchase Order Number  
Purchase Order Date  
QMCIL Reference  
Serial Number

### On receipt of your shipment

Please check that your equipment has arrived safely. Please advise QMC Instruments if you suspect any damage has been incurred during transport and delivery or if any of the items are missing.

This operating manual contains instructions for operation of the cryostat, together with QMC Instruments Ltd. test performance data, against which our guarantee is given as stated above. The user is advised to read this document carefully prior to operation of the cryostat and is reminded that our guarantee will be invalidated if the equipment is damaged through misuse.

Signed.....  
Ken Wood - Director, QMC Instruments Ltd.

Date.....

QMC Instruments technical staff will be happy to advise you if you have any questions or difficulties. The contact details are:

Ken Wood (Sales and Marketing Director)  
QMC Instruments Ltd  
Cardiff University  
School of Physics and Astronomy  
Queens Buildings  
The Parade  
Cardiff  
CF24 3AA  
UK

T. +44 (0) 29 2045 1071  
F. +44 (0) 29 2045 1271  
E. k.wood@terahertz.co.uk  
www.terahertz.co.uk